
4/1990

Gebrauchsanweisung
Instruction Sheet
Mode d'emploi

516 20

516 20 mit 300 01
516 20 with 300 01
516 20 avec 300 01

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Accessories for demonstrating Cou-
lomb's law of electrostatics, and
forces set up in an electric field

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1. Description

The following parts are supplied under Catalogue Number 516 20 as accessories for use with the torsion balance (516 01) to measure electrical forces:

- 1 Metal sphere, dia. 30 mm, mounted on a highly insulating holder, for insertion into the system of the torsion balance.
- 2 Metal spheres, dia. 30 mm, mounted on highly insulating rods approx. 21 cm in length.
- 1 Handle, approx. 33 cm long, with hole for insertion of a metal sphere with rod.
- 1 Adjustable stand with millimetre scale with hole for insertion of a metal sphere with rod. The stand can be adjusted to spacings ranging anywhere between 0 and 25 cm. It is provided with a stand rod of approx. 28 cm length for mounting in a stand base (e. g. 300 01).
- 1 Electrostatic probe plate, an aluminium disc of 40 mm dia. and 0.5 mm thickness with highly insulating holder, for insertion into the system of the torsion balance.

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2. Experiments

1. Coulomb's law of electrostatics

see New Physics Leaflets for Colleges and Universities, Vol. 1 (599 952)

2. Force exerted on a charge in a homogeneous electric field ($F = E Q$)

The electrostatic probe plate is inserted into the system of the torsion balance and placed at the centre of the parallel plate capacitor (544 22) at maximum plate spacing. As a source of potential for the field of the parallel plate capacitor, and to charge the electrostatic probe plate, the high tension power supply (522 37) is suitable.

When performing this experiment, the following points should be noted:

The mid-point of the high tension supply to the parallel plate capacitor must be earthed.

The electrostatic probe plate must be in the exact centre of the parallel plate capacitor. (If the plates are earthed and the probe plate is charged, there must be no deflection.)

While charging the probe plate, both capacitor plates must be earthed, or there will be interferences due to electrostatic induction.

The potential difference across the parallel plate capacitor must not be too small (not less than 2 kV) and the voltage for charging the probe plate must not be too large (not more than 1 kV) or the electrostatic image force will interfere with the determinations.

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