5/1990 -A-

Gebrauchsanweisung  
Instruction Sheet  
Mode d'emploi

546 00  
546 50  
546 02


Fig. 1

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The Wulf electrostatic projector is an apparatus which is used for many purposes, ranging from the simplest electrostatic experiments to measurements on radioactive substances. Each experiment can be rendered visible to a large circle of students by means of the projection attachment or an experiment projector (e. g. FANTEX multi-purpose projector, 443 10), or optical accessories.

## 1 Description

The case of the Wulf electrostatic projector (546 00) is sealed at the front and at the rear by glass plates, and has two openings, one at the top and the other at one side. An adjustable counter-electrode (E) (Fig. 2), is fitted in the lateral opening and insulated against the case. The electrode is held in its respective position by the screw (A), and the adjustment is rendered reproducible by means of the setting ring (B). For cleaning, the insulator and counter-electrode can be removed after slackening 2 grub screws.

A large insulator which can be taken out is fitted in the top opening, and this is provided with a boring to carry the system (S). The bottom of the case is provided with a 4 mm boring for earthing ( $K_4$ ) and two grooves for locating the electrostatic projector in the rods of an experiment projector.

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The bottom of the case is also provided with a shaft, 10 mm  $\phi$  and 10 cm long, which can be screwed off so that the apparatus can be fastened in an extension case supplied with it and secured with the latter in a stand base (300 02) or in a boss. The projection attachment can be fitted on this shaft, too. The shaft is provided with a 4 mm boring for earthing ( $K_3$ ).

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The system is assembled as follows (Fig. 3):

One end of a thin aluminium leaf is secured directly to the system carrier, the other via a fine quartz loop. It can be seen extremely well when projected, since it is turned at an angle of  $180^\circ$  between its ends.

As a result of its special arrangement, which can be seen in Fig. 3, the loop allows the centre part of the leaf to deflect up to 20 mm. The deflection for a certain charge depends upon the position of the counter-electrode. The sensitivity of the Wulf electrostatic increases somewhat quadratically to the laid-on voltage until a labile position is finally reached. With each further increase in voltage, the leaf strikes the counter-electrode and discharge.

These characteristics of its system render the Wulf electrostatic suitable for special measuring methods (see Section 3). Other advantages of the system are its almost inertialess adjustment, its reproducibility, and its insensitiveness to vibrations.

A spare system in an aluminium box is supplied with each electrostatic. Additional spare systems can be obtained under Cat. No. 546 50. A capacitor plate with a 4 mm plug is also supplied with the Wulf electrostatic. When carrying out electrostatic experiments, this plug should be inserted in socket  $K_1$ .


Fig. 3

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
Fig. 4

The projection attachment (546 02) is used for projection if no projector is being used, and consists of two short stand rods located axially to each other and held by a ring-shaped boss with socket, and an adjusting screw. A projecting lens is slid over the one column, and an adjustable cross clamp over the other to secure the lamp housing.

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## 2 Assembly for projection

### a) Projection by means of the projection attachment

The projection attachment (546 02) is slid on the shaft of the electroscope, and the electroscope is secured in a stand base (300 02) either directly or via a stand rod, 25 cm long (300 41), by means of a boss or clamp (301 01 or 03) (Fig. 4). To illuminate the apparatus, place the lamp housing (450 60) with fitted condenser with diaphragm holder (460 20) in the cross-clamp and slide it right up to the electroscope. The lens, electroscope and lamp should all be located at the same level; the projection screen can be placed 2 to 8 metres away. The well defined image of the system is produced by sliding the lens. Finally the lamp filament should be placed in the vertical, and the lamp should be moved in the housing until the lens is uniformly illuminated.

### b) Projection by means of an experiment projector (Fig. 5)

The shaft of the electroscope will be screwed off, and the electroscope will be secured to the rods of the projector (e.g. FANTEX multi-purpose projector, 443 10)

### c) Projection by means of optical accessories (Fig. 6)

The Wulf electroscope is assembled with the aid of 1 stand base (300 02), 2 rods, 25 cm long (300 41), and 4 Leybold bosses (301 01), as shown in Fig. 6.

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The lamp-house with lamp (450 60/51) and the condenser (460 20) are suited for illumination; for projection, a lens  $f=100$  mm or  $f=200$  mm (460 03 or 04) can be used, dependent upon the distance of the projection wall. This assembly can also be made on an optical bench (e.g. 460 43).

### 3 Operation

The system is in the correct position when the space between the carrier and the leaf can be clearly seen.

The Wulf electrostatic is connected to the mains by means of three sockets (see Fig. 2):  $K_1$  is connected to the system,  $K_2$  to the counter-electrode, and  $K_3$  or  $K_4$ , respectively, to the case.

The Wulf electrostatic is used for three measuring methods:

#### a) Circuit without auxiliary charge

The counter-electrode is earthed by connecting  $K_2$  to  $K_3$  or  $K_4$ , respectively. The voltage which is to be observed (e. g. that of the sliding parallel-plate capacitor 544 22) is laid on between the system ( $K_1$ ) and the earth. The system can also be charged either by contact or by influence.

When the charge reaches a certain value, dependent upon the position of the counter-electrode, the lead of the system will deflect against the counter-electrode, and will be discharged again. The deflection of the leaf can also be utilised for charging by influence. A rubbed rod made of glass, hard rubber or plastic material should be placed close up until the leaf deflects. The counter-electrode should then be withdrawn a little, and the rod removed.

With the circuit without auxiliary charge the measuring range is between 200 and 1000 V, dependent upon the position of the counter-electrode. When projecting, and when the counter-electrode is approached to the system, exact measurements of up to 20 V can be made.

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#### b) Circuit with auxiliary charge

The sensitivity of the leaf always reaches a maximum just prior to the leaf deflecting against the counter-electrode. Therefore, the loaded, insulated system is also suited for measuring small voltages of about 1 V. The system ( $K_1$ ) should be loaded either by contact or by influence, and the voltage which is to be measured should be laid on either between ( $K_2$ ) and ( $K_3$ ) or ( $K_4$ ), respectively, or with the counter-electrode charged between ( $K_1$ ) and ( $K_3$ ) or ( $K_4$ ), respectively. With this circuit, it is possible, for example, to demonstrate an induced voltage pulse in a coil electrostatically.

#### c) Drop method

This method is not designed to measure a fixed deflection but to observe the leaf striking against the counter-electrode. Since the charge which is taken away is the same for each deflection, the number of deflections in a certain space of time gives the measurement of the amperage. The method is particularly suited for demonstrating the ionization of the air by radioactive rays, by the photoeffect, or by flame gas. The amperage is a measurement for the number of ions formed.

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### 4 Precautionary measures

The system of the Wulf electroscopie is almost insensitive to vibrations, so that it is not necessary to take any special precautions when transporting the apparatus. On the other hand, we would draw your attention to several points which experience has shown to be the cause of damage to the system.

- a) When taking the system out of the Wulf electroscopie, the system may be easily damaged due to the counter-electrode being close to it. These manipulations should, therefore, only be carried out with special care and when the counter-electrode has been withdrawn.
- b) Also when mounting and removing the supplementary apparatus, the system may be drawn out of the case and be damaged. In order to avoid this, the system should always be held with one hand by its insulator.
- c) When bringing the counter-electrode closer, care should be taken that the system is in the correct position.

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- d) When carrying out electrostatic experiments in which the leaf deflects frequently, the counter-electrode must be earthed, otherwise it will become charged with the result that the newly charged leaf will be drawn to the case and not to the counter-electrode, thus being exposed to damage.
- e) In all experiments where a productive source of current is used (power unit, capacitor), a safety resistor of at least 100 k $\Omega$  must at all costs be connected up in series. Well suited for this purpose is the safety resistor 100 M $\Omega$  (536 22).

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## 5 Installation of a spare system

The installation of a spare system is by no means a difficult process. Withdraw the counter-electrode, slightly turn the insulator with the system and remove it from the case, unscrew the top pin (Fig. 7), and remove the support of the damaged system. Only after having completed these preparations, take the spare system out of the box, and then screw it tightly into the boring of the insulator. When placing it into position, take care not to knock it against the casing because the protruding quartz loop is very sensitive. The system should be projected and, as described above, located in the correct position. The support of the damaged system should be screwed tightly in the case provided, and can be sent either to our works at Hürth or to one of our agents to have a new leaf fitted.

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